

TECH SPECS

- ✓ **MODEL** – VITALS360® Device VC-001
- ✓ **AC ADAPTER** – LXCP12-005200AGF
- ✓ **TYPE C CABLE** – VC-0017
- ✓ **CPU** – Qualcomm® MSM8909
- ✓ **GPU** – Adreno® 306
- ✓ **MEMORY** – 1GB RAM LPDDR4 / 8GB ROM LPDDR4
- ✓ **OPERATING SYSTEM** – Android®
- ✓ **SIZE** – Approximately L 5.71inch x W 3.15inch x T 0.98inch (145mm x 80mm x 25mm)
- ✓ **WEIGHT** – Gross Weight of Complete package 0.552lb (250grams) VITALS360® Device 0.155lb (70grams) Arm Cuff 0.143lb (65grams)
- ✓ **SPLASH, WATER, AND IP22** – Under IEC 60529
- ✓ **DUST RESISTANT CONNECTOR** – Type-C connector support
- ✓ **DISPLAY** – 3.66inch (diagonal) all-screen LED
- ✓ **MULTI-TOUCH DISPLAY** – 800-by-480-pixel resolution at 458 ppi
2,000,000:1 contrast ratio (typical)
- ✓ **CAMERA** – Rear Camera: 5MP PDAF. Front Camera: 2MP FF
- ✓ **MIC** – Yes
- ✓ **SPEAKER** – Yes , 55dBA – 80dBA at a Distance of 10cm
- ✓ **BATTERY TYPE** – Non-removable Size: 1500mAh
- ✓ **CHARGING TIME** – 5 Hours
- ✓ **BATTERY RECYCLE LIFE** – >300times
- ✓ **POWER SOURCE** – 3.7V 1500 mAh or AC adapter (optional, INPUT AC100-240V 50/60Hz
0.5A max)
- ✓ **LANGUAGES** – English (Australia, UK, U.S.), Chinese (Simplified), French (Canada, France),
Spanish (Latin America, Mexico, Spain), Arabic
- ✓ **KEYBOARD SUPPORT** – English (Australia, UK, U.S.), Chinese (Simplified), French (Canada,
France), Spanish (Latin America, Mexico, Spain), Arabic



vitals360
Improving Efficiency

TECH SPECS

