

VoCare offers over the air (OTA) Firmware and Application layer updates via our Mobile Device Management partner Miradore. The Vitals360® device also has the ability to update device settings and bulk deployment option policy for enterprise customers. VoCare MDM Reduces deployment time and ensure the consistency of your end-user environments. Maximizes uptime with increase efficiency and ensures consistency across all devices.

VoCare Mobile Device Management Solutions

VoCare Mobile Device Management

The VoCare Vitals360® deployments are pre-configured to your organization specifications and allows for client configuration change management. VoCare Mobile Device Management solution allows users to update their configuration for use in remote facilities or locations over the air (OTA) as needed, by group or individually.

The VoCare Health status application sends crucial analytics and device health statistics to easily identify and mitigate any operational issues that may arise over time; such as battery degradation or cellular and Wi-Fi connectivity issues. This device health feature enables VoCare Vitals360® to have a long issue free life-span.

VoCare (MDM) for system configuration:

Consistency is crucial to effectively support end users. We can help ensure your systems have the correct settings and versions of applications. To leverage the Mobile Device Management installed on the device, VoCare creates a Business Compliance Plan to synchronize your devices concurrently and automatically, to reduce deployment time.

This plan ensures every employee or customer Vitals360® looks proper to your IT staff. These management tools make it easier to achieve a consistent, efficient environment.

Mobile Device Management Pre-Configuration.

VoCare MDM delivers Pre-configured, ready to use Vitals360® devices faster and easier. Eliminate hours of inefficient end user and IT support time to get the Vitals360® device set up, connected and operational. Utilize the VoCare MDM solution to quickly and accurately deploy hundreds to thousands of devices with identical configuration, and company required modifications.

The Mobile Device Management Connected Provisioning helps provide a custom configured Vitals360® device to an end user faster with less interaction. VoCare leverages the cloud to speed the way devices are configured and managed. This company custom solution is designed to benefit IT policy and the end user application

